

Test Definition: FUNBL

Fungitell, BAL

Overview

Method Name

Limulus Amebocyte Lysate (LAL) Pathway

NY State Available

Yes

Specimen

Specimen Type

Lavage

Specimen Required

Specimen Type: BAL

Sources: Bronchoalveolar lavage **Container/Tube**: Sterile container **Specimen Volume**: 1 – 3 mL

Note: Ship frozen

Specimen Minimum Volume

0.5 mL

Reject Due To

Hemolysis	Mild reject; Gross reject	
Lipemia	Mild reject; Gross reject	
Icterus	Mild reject; Gross reject	
Other	NA	

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Lavage	Ambient	4 days	
	Refrigerated	7 days	
	Frozen (preferred)	30 days	

Clinical & Interpretive



Test Definition: FUNBL

Fungitell, BAL

Clinical Information

The Fungitell Beta – D Glucan assay detects (1,3) – Beta-D-glucan from the following pathogens: Candida spp., Acremonium, Aspergillus spp., Coccidioides immitis, Fusarium spp., Histoplasma capsulatum, Trichosporon spp., Sporothrix schenckii, Saccharomyces cerevisiae, and Pneumocystis jiroveci.

The Fungitell Beta-D Glucan assay does not detect certain fungal species such as the genus Cryptococcus, which produces very low levels of (1,3) $\hat{a} \in Beta-D-glucan$, nor the Zygomycetes, such as Absidia, Mucor, and Rhizopus, which are not known to produce (1,3) $\hat{a} \in Beta-D-glucan$. Studies indicate Blastomyces dermatitidis is usually not detected due to little (1,3) $\hat{a} \in Beta-D-glucan$ produced in the yeast phase.

Reference Values

A reference range for specimens other than serum has not been established.

Interpretation

The performance characteristics of the Fungitell assay in BAL have been determined by Viracor-Eurofins; there are no established criteria for the interpretation of Fungitell results from BAL fluid. Research studies have evaluated the use of the Fungitell assay in BAL in both immunocompromised patients (Mycopathologia (2013) 175:33-41) and acute eosinophilic pneumonia (Chest (2003) 123:1302-1307).

Performance

PDF Report

No

Day(s) Performed

Monday through Saturday

Report Available

1 to 6 days

Performing Laboratory Location

Eurofins Viracor

Fees & Codes

Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

CPT Code Information

87449



Test Definition: FUNBL

Fungitell, BAL

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FUNBL	Fungitell, BAL	Not Provided

Result ID	Test Result Name	Result LOINC® Value
FUNBL	Fungitell, BAL	Not Provided