

Overview

Useful For

Detecting drug use involving stimulants, benzodiazepines, and opioids

This test is **not intended** for use in employment-related testing.

Profile Information

| Test Id | Reporting Name | Available Separately | Always Performed |
|---------|------------------------------------|----------------------|------------------|
| LPCM | List Patient's Current Medications | No | Yes |
| ADULT | Adulterants Survey, U | Yes | Yes |
| TOPSU | Targeted Opioid Screen, U | Yes, (order TOSU) | Yes |
| TABSU | Targeted Benzodiazepine Screen, U | Yes, (order TBSU) | Yes |
| TSTIM | Targeted Stimulant Screen, U | Yes, (order TSPU) | Yes |

Testing Algorithm

Testing begins with an adulterant survey. If the sample is found to be adulterated, testing will end, and the remaining tests will be canceled.

If the specimen is normal or only diluted, testing will proceed.

Method Name

ADULT: Spectrophotometry

TOPSU, TABSU, TSTIM: Liquid Chromatography Tandem Mass Spectrometry, High-Resolution Accurate Mass (LC-MS/MS HRAM)

NY State Available

Yes

Specimen

Specimen Type

Urine

Ordering Guidance

This test does not screen for drug classes other than those listed in Reference Values.

Specimen Required

Supplies: Urine Tubes, 10 mL (T068)

Collection Container/Tube: Plastic urine container

Submission Container/Tube: Plastic, 10 mL tube

Specimen Volume: 5 mL

Collection Instructions:

1. Collect a random urine specimen.
2. Submit 5 mL in 1 plastic bottle.
3. No preservative

Additional Information:

1. No specimen substitutions.
2. Submitting less than 5 mL may compromise the ability to perform all necessary testing.
3. STAT requests are **not accepted** for this test.

Forms

If not ordering electronically, complete, print, and send a [Therapeutics Test Request](#) (T831) with the specimen.

Specimen Minimum Volume

2 mL

Reject Due To

All specimens will be evaluated at Mayo Clinic Laboratories for test suitability.

Specimen Stability Information

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------------|---------|-------------------|
| Urine | Refrigerated (preferred) | 14 days | |
| | Frozen | 14 days | |

Clinical & Interpretive**Clinical Information**

The targeted opioid, benzodiazepine, and stimulant screen portions are performed by liquid chromatography tandem mass spectrometry, high-resolution accurate mass (LC-MS/MS HRAM) and are completed for all opioids, benzodiazepines, and stimulants.

Opioids are a large class of medications commonly used to relieve acute and chronic pain or help manage opioid abuse and dependence. Medications that fall into this class include buprenorphine, codeine, fentanyl, hydrocodone, hydromorphone, methadone, morphine, oxycodone, oxymorphone, tapentadol, tramadol, and others. Opioids work by

binding to the opioid receptors that are found in the brain, spinal cord, gastrointestinal tract, and other organs. Common side effects of opioids include drowsiness, confusion, nausea, constipation, and, in severe cases, respiratory depression. These are dose dependant and vary with tolerance. These medications can also produce physical and psychological dependence and have a high risk for abuse and diversion, which is one of the main reasons many professional practice guidelines recommend compliance testing in patients prescribed these medications.

Opioids are readily absorbed from the gastrointestinal tract, nasal mucosa, lungs, and after subcutaneous or intramuscular injection. Opioids are primarily excreted from the kidney in both free and conjugated forms. This assay does not hydrolyze the urine sample and looks for both parent drugs and metabolites (including glucuronide forms). The detection window for most opioids in urine is approximately 1 to 3 days with longer detection times for some compounds (eg, methadone).

Benzodiazepines represent a large family of medications used to treat a wide range of disorders from anxiety to seizures and are also used in pain management. With a high risk for abuse and diversion, professional practice guidelines recommend compliance monitoring for these medications using urine drug tests. However, traditional benzodiazepine immunoassays suffer from a lack of cross-reactivity with all the benzodiazepines, so many compliant patients taking either clonazepam (Klonopin) or lorazepam (Ativan) may screen negative by immunoassay but are positive when confirmatory testing is done. The new targeted benzodiazepine screening test provides a more sensitive and specific test to check for compliance to all the commonly prescribed benzodiazepines and looks for both parent drug and metabolites in the urine.

Stimulants are sympathomimetic amines that stimulate the central nervous system activity and, in part, suppress the appetite. Amphetamine and methamphetamine are also prescription drugs used in the treatment of narcolepsy and attention-deficit disorder/attention-deficit hyperactivity disorder (ADHD). Methylphenidate is another stimulant used to treat ADHD. Phentermine is indicated for the management of obesity. All other amphetamines (eg, methylenedioxymethamphetamine: MDMA) are Drug Enforcement Administration-scheduled Class I compounds. Due to their stimulant effects, the drugs are commonly sold illicitly and abused. Physiological symptoms associated with very high amounts of ingested amphetamine or methamphetamine include elevated blood pressure, dilated pupils, hyperthermia, convulsions, and acute amphetamine psychosis.

This test is intended to be used in a setting where the test results can be used to make a definitive diagnosis.

Reference Values**ADULTERANT SURVEY:**

Cutoff concentrations

Oxidants: 200 mg/L

Nitrites: 500 mg/L

TARGETED OPIOID SCREEN:

Not detected (Positive results are reported with qualitative "Present" results)

Cutoff concentrations:

Codeine: 25 ng/mL

Codeine-6-beta-glucuronide: 100 ng/mL
Morphine: 25 ng/mL
Morphine-6-beta-glucuronide: 100 ng/mL
6-Monoacetylmorphine: 25 ng/mL
Hydrocodone: 25 ng/mL
Norhydrocodone: 25 ng/mL
Dihydrocodeine: 25 ng/mL
Hydromorphone: 25 ng/mL
Hydromorphone-3-beta-glucuronide: 100 ng/mL
Oxycodone: 25 ng/mL
Noroxycodone: 25 ng/mL
Oxymorphone: 25 ng/mL
Oxymorphone-3-beta-glucuronide: 100 ng/mL
Noroxymorphone: 25 ng/mL
Fentanyl: 2 ng/mL
Norfentanyl: 2 ng/mL
Meperidine: 25 ng/mL
Normeperidine: 25 ng/mL
Naloxone: 25 ng/mL
Naloxone-3-beta-glucuronide: 100 ng/mL
Methadone: 25 ng/mL
2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP): 25 ng/mL
Propoxyphene: 25 ng/mL
Norpropoxyphene: 25 ng/mL
Tramadol: 25 ng/mL
O-desmethyltramadol: 25 ng/mL
Tapentadol: 25 ng/mL
N-desmethiltapentadol: 50 ng/mL
Tapentadol-beta-glucuronide: 100 ng/mL
Buprenorphine: 5 ng/mL
Norbuprenorphine: 5 ng/mL
Norbuprenorphine glucuronide: 20 ng/mL

TARGETED BENZODIAZEPINE SCREEN:

Not detected (Positive results are reported with qualitative "Present" results)

Cutoff concentrations:

Alprazolam: 10 ng/mL
Alpha-hydroxyalprazolam: 10 ng/mL
Alpha-hydroxyalprazolam glucuronide: 50 ng/mL
Chlordiazepoxide: 10 ng/mL
Clobazam: 10 ng/mL
N-desmethyloclobazam: 200 ng/mL

Clonazepam: 10 ng/mL
7-Aminoclonazepam: 10 ng/mL
Diazepam: 10 ng/mL
Nordiazepam: 10 ng/mL
Flunitrazepam: 10 ng/mL
7-Aminoflunitrazepam: 10 ng/mL
Flurazepam: 10 ng/mL
2-Hydroxy ethyl flurazepam: 10 ng/mL
Lorazepam: 10 ng/mL
Lorazepam glucuronide: 50 ng/mL
Midazolam: 10 ng/mL
Alpha-hydroxymidazolam: 10 ng/mL
Oxazepam: 10 ng/mL
Oxazepam glucuronide: 50 ng/mL
Prazepam: 10 ng/mL
Temazepam: 10 ng/mL
Temazepam glucuronide: 50 ng/mL
Triazolam: 10 ng/mL
Alpha-hydroxytriazolam: 10 ng/mL
Zolpidem: 10 ng/mL
Zolpidem phenyl-4-carboxylic acid: 10 ng/mL

TARGETED STIMULANT SCREEN:

Not detected (Positive results are reported with qualitative "Present" results)

Cutoff concentrations:

Methamphetamine: 100 ng/mL
Amphetamine: 100 ng/mL
3,4-Methylenedioxymethamphetamine (MDMA): 100 ng/mL
3,4-Methylenedioxy-N-ethylamphetamine (MDEA): 100 ng/mL
3,4-Methylenedioxyamphetamine (MDA): 100 ng/mL
Ephedrine: 100 ng/mL
Pseudoephedrine: 100 ng/mL
Phentermine: 100 ng/mL
Phencyclidine (PCP): 20 ng/mL
Methylphenidate: 20 ng/mL
Ritalinic acid: 100 ng/mL

Interpretation

A positive result derived by this testing indicates that the patient has used one of the drugs detected by these techniques in the recent past.

For information about drug testing, including estimated detection times and [Result Interpretations](#), see [Controlled](#)

[Substance Monitoring](#) on MayoClinicLabs.com.

Cautions

No significant cautionary statements

Clinical Reference

1. Physicians' Desk Reference: 60th ed. Medical Economics Company; 2006
2. Brunzman LL Lazo JS, Parker KL, eds. Goodman and Gilman's: The Pharmacological Basis of Therapeutics. 11th ed. McGraw-Hill Book Company; 2006
3. Langman LJ, Bechtel LK, Holstege CP. Clinical toxicology. In: Rifai N, Chiu RWK, Young I, Burnham CAD, Wittwer CT, eds. Tietz Textbook of Laboratory Medicine. 7th ed. Elsevier; 2023:chap 43
4. Gutstein HB, Akil H: Opioid analgesics. In: Brunton LL, Lazo JS, Parker KL, eds. Goodman and Gilman's: The Pharmacological Basis of Therapeutics. 11th ed. McGraw-Hill Companies; 2006:chap 21
5. Rovine T, Ferrero CL, American Pain Society: Chronic Pain in America: Roadblocks to Relief. Roper Starch Worldwide, Inc; 1999. Updated 2001. Accessed December 12, 2024. Available at <http://accurateclinic.com/wp-content/uploads/2016/04/Chronic-Pain-In-America-Roadblocks-To-Relief-1999.pdf>
6. Magnani B, Kwong T. Urine drug testing for pain management. Clin Lab Med. 2012;32(3):379-390
7. Jannetto PJ, Bratanow NC, Clark WA, et al. Executive summary: American Association of Clinical Chemistry Laboratory Medicine Practice Guideline-using clinical laboratory tests to monitor drug therapy in pain management patients. J Appl Lab Med. 2018;2(4):489-526
8. McMillin GA, Marin SJ, Johnson-Davis KL, Lawlor BG, Strathmann FG. A hybrid approach to urine drug testing using high-resolution mass spectrometry and select immunoassays. Am J Clin Pathol. 2015;143(2):234-240
9. Cone EJ, Caplan YH, Black DL, Robert T, Moser F. Urine drug testing of chronic pain patients: licit and illicit drug patterns. J Anal Toxicol. 2008;32(8):530-543

Performance**Method Description**

Adulterant:

All results are measured using spectrophotometry at wavelengths specified by the reagent manufacturer. The use of a refractometer may also be used in the specific gravity measurement.(Package inserts: Specimen Validity Test Creatinine. Roche Diagnostics; V3.0, 08/2015; Specimen Validity Test Nitrite. Roche Diagnostics; V3.0, 08/2018, Specimen Validity Test Oxidant. Roche Diagnostics; V 3.0, 08/2018; Specimen Validity Test pH Roche Diagnostics; V3.0, 02/2019, Specimen Validity Test Specific Gravity. Roche Diagnostics; V4.0, 08/2022)

Targeted Screening Panels for opioids, benzodiazepines, and stimulants:

The urine sample is diluted with internal standard and clinical laboratory reagent water and then analyzed by liquid chromatography tandem mass spectrometry using a high-resolution accurate mass orbitrap detector.(Unpublished Mayo method)

PDF Report

Test Definition: CSMTU

Controlled Substance Monitoring Targeted Profile, 17 Drug Classes, Mass Spectrometry, Random, Urine

No

Day(s) Performed

Monday through Saturday

Report Available

3 to 4 days

Performing Laboratory Location

Rochester

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. It has not been cleared or approved by the US Food and Drug Administration.

CPT Code Information

80347
80364
80326
G0482 (if appropriate)

LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---------|-------------------------------------|--------------------|
| CSMTU | CSM Targeted Drug Profile,17,HRMS,U | 69739-1 |

| Result ID | Test Result Name | Result LOINC® Value |
|-----------|------------------|---------------------|
| 20606 | Creatinine, U | 2161-8 |
| 22312 | Specific Gravity | In Process |
| 23509 | pH | 2756-5 |
| 23511 | Oxidants | 58714-7 |
| 23510 | Nitrites | 32710-6 |
| 30914 | Comment | 48767-8 |
| 42323 | Codeine | 19411-8 |

Test Definition: CSMTU

Controlled Substance Monitoring Targeted
Profile, 17 Drug Classes, Mass Spectrometry,
Random, Urine

| | | |
|--------|--|---------|
| 42324 | Codeine-6-beta-glucuronide | 89310-7 |
| 42325 | Morphine | 19597-4 |
| 42326 | Morphine-6-beta-glucuronide | 89308-1 |
| 42327 | 6-monoacetylmorphine | 19321-9 |
| 42328 | Hydrocodone | 19482-9 |
| 42329 | Norhydrocodone | 89304-0 |
| 42330 | Dihydrocodeine | 19446-4 |
| 42331 | Hydromorphone | 19486-0 |
| 42332 | Hydromorphone-3-beta-glucuronide | 89309-9 |
| 42333 | Oxycodone | 19642-8 |
| 42334 | Noroxycodone | 89303-2 |
| 42335 | Oxymorphone | 19646-9 |
| 42336 | Oxymorphone-3-beta-glucuronide | 89301-6 |
| 42337 | Noroxymorphone | 89302-4 |
| 42338 | Fentanyl | 59673-4 |
| 42339 | Norfentanyl | 43199-9 |
| 42340 | Meperidine | 19532-1 |
| 42341 | Normeperidine | 27920-8 |
| 42342 | Naloxone | 42618-9 |
| 42343 | Naloxone-3-beta-glucuronide | 89307-3 |
| 42344 | Methadone | 19550-3 |
| 42345 | EDDP | 93495-0 |
| 42346 | Propoxyphene | 19429-0 |
| 42347 | Norpropoxyphene | 19632-9 |
| 42348 | Tramadol | 19710-3 |
| 42349 | O-desmethyltramadol | 86453-8 |
| 42350 | Tapentadol | 72485-6 |
| 42351 | N-desmethyltapentadol | 89306-5 |
| 42352 | Tapentadol-beta-glucuronide | 89300-8 |
| 42353 | Buprenorphine | 93494-3 |
| 42354 | Norbuprenorphine | 82371-6 |
| 42355 | Norbuprenorphine glucuronide | 89305-7 |
| 65059 | Opioid Interpretation | 69050-3 |
| 604871 | Alprazolam | 94116-1 |
| 604867 | Alpha-Hydroxyalprazolam | 19325-0 |
| 604891 | Alpha-Hydroxyalprazolam Glucuronide | 94115-3 |
| 604872 | Chlordiazepoxide | 19385-4 |
| 604889 | Clobazam | 94114-6 |
| 604890 | N-Desmethyloclobazam | 94113-8 |
| 604873 | Clonazepam | 19399-5 |

Test Definition: CSMTU

Controlled Substance Monitoring Targeted
Profile, 17 Drug Classes, Mass Spectrometry,
Random, Urine

| | | |
|--------|--|---------|
| 604267 | 7-aminoclonazepam | 94112-0 |
| 604874 | Diazepam | 19443-1 |
| 604880 | Nordiazepam | 19624-6 |
| 604875 | Flunitrazepam | 19466-2 |
| 604866 | 7-aminoflunitrazepam | 94111-2 |
| 604876 | Flurazepam | 19474-6 |
| 604868 | 2-Hydroxy Ethyl Flurazepam | 94110-4 |
| 604877 | Lorazepam | 19520-6 |
| 604878 | Lorazepam Glucuronide | 94109-6 |
| 604879 | Midazolam | 19585-9 |
| 604869 | Alpha-Hydroxy Midazolam | 94108-8 |
| 604881 | Oxazepam | 19638-6 |
| 604882 | Oxazepam Glucuronide | 94107-0 |
| 604883 | Prazepam | 19678-2 |
| 604884 | Temazepam | 19698-0 |
| 604885 | Temazepam Glucuronide | 94106-2 |
| 604886 | Triazolam | 19714-5 |
| 604870 | Alpha-Hydroxy Triazolam | 94105-4 |
| 604887 | Zolpidem | 94104-7 |
| 604888 | Zolpidem Phenyl-4-Carboxylic acid | 94103-9 |
| 604949 | Benzodiazepine Interpretation | 69050-3 |
| LPCM | List Patient's Current Medications | 66423-5 |
| 610273 | Methamphetamine | 19554-5 |
| 610274 | Amphetamine | 19343-3 |
| 610275 | 3,4-methylenedioxyamphetamine (MDMA) | 19568-5 |
| 610276 | 3,4-methylenedioxy-N-ethylamphetamine (MDEA) | 59844-1 |
| 610277 | 3,4-methylenedioxyamphetamine (MDA) | 19565-1 |
| 610278 | Ephedrine | 99108-3 |
| 610279 | Pseudoephedrine | 99109-1 |
| 610280 | Phentermine | 19674-1 |
| 610281 | Phencyclidine (PCP) | 19659-2 |
| 610282 | Methylphenidate | 19577-6 |
| 610283 | Ritalinic acid | 99110-9 |
| 610284 | Stimulant Interpretation | 54247-2 |