

## Overview

**Method Name**

Direct EIA/ELISA

**NY State Available**

No

## Specimen

**Specimen Type**

Urine

**Specimen Required**

**Patient Preparation:** Patient should not be on aspirin, indomethacin, or anti-inflammatory medications, if possible, for at least 48 hours prior to collection of specimen. Specimen must be frozen within 30 minutes of collection.

**Specimen Type:** Urine

**Submission Container/Tube:** Plastic, 10-mL tube (T068)

**Specimen Volume:** 10 mL

**Collection Instructions:**

1. Collect random urine (NO preservative).
3. Freeze immediately and send specimen frozen in the plastic, 10-mL urine tube (T068)

**Note:** 24 hours urine collection is not acceptable.

**Specimen Minimum Volume**

5 mL

**Reject Due To**

Hemolysis	NA
Lipemia	NA
Icteric	NA
Other	NA

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Urine	Frozen	180 days	

## Clinical & Interpretive

### Clinical Information

Prostaglandins are fatty acids derived from arachidonic acid metabolism. They are closely related to the Thromboxanes and Leukotrienes. Prostaglandin D2 is derived mainly from Prostaglandin H2 and is metabolized to Dihydroketo Prostaglandin D2. Prostaglandin D2 is excreted directly into the urine. The sites of highest Prostaglandin D2 activity are the brain, spinal cord, intestines, and stomach. Prostaglandin D2 is the major Prostaglandin produced by uterine tissue. Prostaglandin D2 is a potent bronchoconstrictor, neuromodulator, and anti-antithrombin agent. It also stimulates the secretion of Pancreatic Glucagon. Prostaglandin D2 has been found to have an anti-metastatic effect on many malignant tumor cells. Prostaglandin D2 production and circulating levels are drastically suppressed by aspirin and indomethacin.

### Reference Values

Up to 175 ng/g Creatinine

## Performance

### PDF Report

Referral

### Day(s) Performed

Upon receipt

### Report Available

10 to 14 days

### Performing Laboratory Location

Inter Science Institute

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### Test Classification

This test was performed using a kit that has not been cleared or approved by the FDA and is designated as research use only. The analytic performance characteristics of this test have been determined by Inter Science Institute. This test is not intended for diagnosis or patient management decisions without confirmation by other medically established means.

---

**CPT Code Information**

84150

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
FPD2U	Prostaglandin D2 (PG D2), U	Not Provided

Result ID	Test Result Name	Result LOINC® Value
FPD2U	Prostaglandin D2 (PG D2), U	Not Provided