

Overview

Method Name

ELISA

NY State Available

Yes

Specimen

Specimen Type

Serum

Specimen Required

Collection Container/Tube: 5 mL Red

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL

Acceptable: SST

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume

0.3 mL (volume does NOT allow for repeat testing)

Reject Due To

Gross hemolysis:	Reject
Thawing:	Warm OK; Cold OK
Gross lipemia:	Reject
Gross icterus	Reject
Other:	Anything other than serum; bacterial contamination

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Ambient	7 days	
	Refrigerated (preferred)	14 days	
	Frozen	60 days	

## Clinical & Interpretive

### Clinical Information

The anti-PM/Scl-100 antibody is associated with younger age, calcinosis and has lower rates of gastrointestinal symptoms, ILD and pulmonary hypertension. There is also evidence of a possibly better survival compared to the presence of either anti-PM/Scl-75 or anti-Scl-70 antibodies.

### Reference Values

Reference Range: <20

### Interpretation:

Negative: <20 units

Weak Positive: 20-39 units

Moderate Positive: 40-80 units

Strong Positive: >80 units

## Performance

### PDF Report

No

### Report Available

16-25 days

### Performing Laboratory Location

Esoterix Endocrinology

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### Test Classification

This test was developed and its performance characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.

### CPT Code Information

86235

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FAPMA	Anti-PM/ScI-100 Ab	101221-0

Result ID	Test Result Name	Result LOINC® Value
FAPMA	Anti-PM/ScI-100 Ab	101221-0