

Overview

Method Name

Real-Time Polymerase Chain Reaction (RT-PCR)

NY State Available

Yes

Specimen

Specimen Type

Plasma

Necessary Information

Specimen source is required.

Specimen Required

Source is required.

Specimen Type: Plasma

Collection Container/Tube: lavender-top (EDTA)

Submission Container Type: 12x75 mm screw capped vial

Specimen Volume: 1 mL

Collection Instructions: Draw blood in a lavender-top (EDTA) tube(s). Spin down and transfer 1 mL EDTA plasma into a plastic, screw capped vial. Send specimen refrigerated.

Specimen Minimum Volume

0.6 mL

Reject Due To

Other	Whole blood or plasma anticoagulated with heparin; calcium alginate swabs, sputum
-------	-----------------------------------------------------------------------------------

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Plasma	Refrigerated (preferred)	7 days	
	Frozen	30 days	
	Ambient	48 hours	

Clinical & Interpretive**Reference Values**

Not Detected

Performance**PDF Report**

No

Day(s) Performed

Monday through Sunday

Report Available

4 to 8 days

Performing Laboratory Location

Quest Diagnostics

Fees & Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its analytical performance characteristics have been determined by Quest Diagnostics. It has not been cleared or approved by FDA. This assay has been validated pursuant to the CLIA regulations and is used for clinical purposes.

CPT Code Information

87799

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FADNA	Adenovirus DNA, Quant RT PCR, P	49340-3

Result ID	Test Result Name	Result LOINC® Value
-----------	------------------	---------------------

Test Definition: FADNA

Adenovirus DNA, Quantitative Real-Time PCR,
Plasma

Z6088	Source	31208-2
Z6089	Adenovirus DNA, QN PCR	49340-3
Z6090	Adenovirus DNA, QN RT PCR	64038-3