

**Overview****Method Name**

High Performance Liquid Chromatography with Ultraviolet Detection (HPLC-UV)

**NY State Available**

Yes

**Specimen****Specimen Type**

Varies

**Specimen Required****Submit only 1 of the following specimens:****Plasma**

Draw blood in a green-top (sodium heparin) tube(s), plasma gel tube is not acceptable. Spin down and send 2 mL sodium heparin plasma refrigerated in a plastic vial.

**Serum**

Draw blood in a plain red-top tube(s), serum gel tube is not acceptable. Spin down and send 2 mL of serum refrigerated in a plastic vial.

**Specimen Minimum Volume**

0.3 mL

**Reject Due To**

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	EDTA

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	14 days	
	Frozen	180 days	

	Ambient	72 hours	
--	---------	----------	--

## Clinical & Interpretive

### Reference Values

Amoxapine	No reference range provided	ng/mL
8-Hydroxyamoxapine	No reference range provided	ng/mL
Combined Total	200 – 400	ng/mL

## Performance

### PDF Report

No

### Day(s) Performed

Monday through Sunday

### Report Available

3 to 7 days

### Specimen Retention Time

2 weeks

### Performing Laboratory Location

Medtox Laboratories, Inc.

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### CPT Code Information

80335

### LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
---------	-----------------	--------------------

## Test Definition: FAMOX

Amoxapine (Asendin) and  
8-Hydroxyamoxapine

FAMOX	Amoxapine (Asendin)	9625-5
-------	---------------------	--------

Result ID	Test Result Name	Result LOINC® Value
Z1120	Amoxapine	3341-5
Z1045	8-Hydroxy-Amoxapine	1698-0
Z1046	Combined Total	9625-5