

Overview**Method Name**

High Performance Liquid Chromatography with Ultraviolet Detection (HPLC-UV)

NY State Available

Yes

Specimen**Specimen Type**

Varies

Specimen Required**Submit only 1 of the following specimens:****Plasma**

Draw blood in a green-top (sodium heparin) tube(s), plasma gel tube is not acceptable. Spin down and send 2 mL sodium heparin plasma refrigerated in a plastic vial.

Serum

Draw blood in a plain red-top tube(s), serum gel tube is not acceptable. Spin down and send 2 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume

0.3 mL

Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	EDTA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Ambient	72 hours	
	Refrigerated (preferred)	14 days	

	Frozen	180 days	
--	--------	----------	--

Clinical & Interpretive

Reference Values

Amoxapine	No reference range provided	ng/mL
8-Hydroxyamoxapine	No reference range provided	ng/mL
Combined Total	200 – 400	ng/mL

Performance

PDF Report

No

Day(s) Performed

Monday through Sunday

Report Available

3 to 7 days

Specimen Retention Time

2 weeks

Performing Laboratory Location

Medtox Laboratories, Inc.

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

80335

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
---------	-----------------	--------------------

Test Definition: FAMOX

Amoxapine (Asendin) and
8-Hydroxyamoxapine

FAMOX	Amoxapine (Asendin)	9625-5
-------	---------------------	--------

Result ID	Test Result Name	Result LOINC® Value
Z1120	Amoxapine	3341-5
Z1045	8-Hydroxy-Amoxapine	1698-0
Z1046	Combined Total	9625-5