

**Overview**

**Method Name**

High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)

**NY State Available**

Yes

**Specimen**

**Specimen Type**

Varies

**Specimen Required**

**\*\*\*Must submit one specimen per order. Specimens cannot be shared between multiple orders.\*\*\***

**Submit only 1 of the following specimens.**

**Serum**

**Specimen Type:** Serum

**Collection Container/Tube:** Red top

**Specimen Volume:** 1 mL

**Collection Instructions:** Draw blood in a plain, red-top tube(s). (Serum gel tube is not acceptable.) Spin down and send 1 mL of serum refrigerated in a plastic preservative-free vial.

**Plasma**

**Specimen Type:** Plasma

**Collection Container/Tube:** Lavender top or pink top (EDTA)

**Specimen volume:** 1 mL

**Collection Instructions:** Draw blood in a lavender-top (EDTA) or Pink top tube(s). (Plasma gel tube is not acceptable.) Spin down and send 1 mL of EDTA plasma refrigerated in a plastic preservative-free vial.

**Specimen Minimum Volume**

0.4 mL

**Reject Due To**

Other	Polymer gel separation tube (SST or PST)
-------	--

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
---------------	-------------	------	-------------------

Varies	Ambient	30 days	
	Refrigerated (preferred)	30 days	
	Frozen	30 days	

## Clinical & Interpretive

### Reference Values

Reporting limit determined each analysis

None Detected ng/mL

Mean steady-state plasma levels following a daily regimen:

80 mg: 26 - 36 ng/mL

160 mg: 52 - 74 ng/mL

320 mg: 154 - 191 ng/mL

## Performance

### PDF Report

No

### Day(s) Performed

Monday through Sunday

### Report Available

8 to 10 days

### Performing Laboratory Location

NMS Labs

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### Test Classification

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the U.S. Food and Drug Administration.

---

**CPT Code Information**

80299

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
FNAD	Nadolol	12406-5

Result ID	Test Result Name	Result LOINC® Value
Z1430	Nadolol	12406-5
Z1867	Reporting Limit	19147-8