



## Eculizumab Monitoring Panel, Serum

**Test ID:** ECMP

**Useful for:**

Monitoring of complement blockage by eculizumab

Assessing the response to eculizumab therapy

Assessing the need for dose escalation

Evaluating the potential for dose de-escalation or discontinuation of therapy in remission states

Monitoring patients who need to be above a certain eculizumab concentration to improve the odds of a clinical response for therapy optimization

This test is **not useful** as the sole basis for a diagnosis or treatment decisions.

**Profile Information:**

Test ID	Reporting Name	Available Separately	Always Performed
EAH50	Eculizumab Complement Blockage, S	No	Yes
ECULI	Eculizumab, S	Yes	Yes
ECUIN	Eculizumab Interpretation, S	No	Yes

**Methods:**

EAH50: Enzyme-Linked Immunosorbent Assay (ELISA)

ECULI: Liquid Chromatography Tandem Mass Spectrometry (LC-MS/MS), High Resolution Accurate Mass

**Reference Values:**

Eculizumab Complement Blockage:

> or =46% normal

Eculizumab:

Lower limit of quantitation =5.0 mcg/mL

>35 Therapeutic concentration for paroxysmal nocturnal hemoglobinuria (PNH)

>50 Therapeutic concentration for atypical hemolytic uremic syndrome (aHUS)

**Specimen Requirements:**

**Patient Preparation:**

1. Fasting preferred.
2. Suggest discontinuing natalizumab at least 4 weeks prior to testing for eculizumab quantitation in serum. Patient should consult the healthcare provider who prescribed this drug to determine if discontinuation is an option. If not, ok to proceed with testing while taking natalizumab.

**Supplies:**

Sarstedt Aliquot Tube, 5 mL (T914)

**Collection Container/Tube:**

Red top (serum gel/SST are **not acceptable**)

**Specimen Volume:**

2 mL in 2 plastic vials, each vial containing 1 mL

**Collection Instructions:**

1. Draw blood immediately before next scheduled dose.
2. Immediately after specimen collection, place the tube on wet ice.
3. After sample has clotted on wet ice, centrifuge at 4 degrees C and aliquot serum into two 5 mL plastic vials.
4. Freeze specimen within 30 minutes of centrifugation. Sample must be placed on dry ice if not frozen immediately.

**Minimum Volume:**

1 mL total in 2 plastic vials, each vial containing 0.5 mL

**Specimen Stability Information:**

Specimen Type	Temperature	Time
Serum	Frozen	14 days
Serum Red	Frozen	14 days

**Cautions:**

Eculizumab complement blockage monitoring is a functional test and is dependent on correct sampling, storage, and shipping conditions. Both degradation by temperature and consumption of complement components will lead to false low function results. These are difficult to differentiate from real complement dysregulation or blockage, and in the event of poor pre-analytical handling, eculizumab concentrations are a more reliable indicator, as not subject to stringent temperature stability.

While pre-analytic handling can lead to falsely low results, it is far less likely that it would lead to falsely normal results.

Complement testing may be ordered in several circumstances where standard treatment includes plasmapheresis or plasma exchange. The procedure itself, if traumatic, may activate complement and therefore, may not be a true reflection of the patient's complement system. The recommendation is to collect blood prior to the plasma exchange whenever possible.

Functional results inconsistent with the clinical history should be verified with a new blood draw.

Specimens should be frozen immediately after collection.

Long term stability is optimal when the sample is kept at -70 degrees Celsius or lower prior to testing.

Results must be interpreted within the clinical context of the patient.

Patients in transition between eculizumab and ravulizumab administration will have a result that is the sum of eculizumab plus ravulizumab in circulation. This assay will not clearly differentiate between these specific analytes and must be interpreted with caution.

Patients actively undergoing therapy with both natalizumab and eculizumab (extremely rare scenario) could present as assay interference. It is suggested patients discuss with their doctors the possibility of discontinuing natalizumab 4 weeks prior to testing. If discontinuation is not possible, it is ok to proceed with testing.

**CPT Code:**

80299

86161

**Day(s) Performed:** Varies

**Report Available:** 3 to 10 days

**Questions**

Contact Amy Ennis, Laboratory Resource Coordinator at 800-533-1710.