

## Adenovirus DNA, Quantitative Real-Time PCR, Plasma

**Test ID:** FADNA

### Explanation:

Due to in-house test availability, test ID: FADNA performed by Quest Diagnostics, will be obsolete effective December 3, 2024.

### Recommended Alternative Test:

## Adenovirus DNA Detection and Quantification, Plasma

**Test ID:** ADVQU

### Methodology:

Real-Time (Quantitative) Polymerase Chain Reaction (rt-qPCR)

### Reference Values:

Undetected

### Specimen Requirements:

<b>Supplies</b>	Sarstedt Aliquot Tube, 5 mL (T914)
<b>Collection Container/Tube:</b>	Lavender top (EDTA)
<b>Submission Container/Tube:</b>	Plastic vial
<b>Specimen Volume:</b>	1.5 mL (plasma)
<b>Collection Instructions:</b>	1. Centrifuge blood collection tube per manufacturer's instructions (eg, centrifuge and aliquot within 2 hours of collection for BD Vacutainer tubes). 2. Aliquot plasma into plastic vial.
<b>Minimum Volume:</b>	0.8 mL

### Specimen Stability Information:

Specimen Type	Temperature	Time
Plasma EDTA	Frozen (preferred)	30 days

	Refrigerated	6 days
--	--------------	--------

**CPT Code:**

87799

**Day(s) Performed:** Tuesday, Friday **Report Available:** 1 to 4 days

**Questions**

Contact MCL Referrals Supervisor Amy Bluhm at 800-533-1710.