

**Reporting Title:** RPR Screen w/ Reflex to Titer,S  
**Performing Location:** Jacksonville

**Ordering Guidance:**  
This test is for monitoring response to therapy in patients treated for syphilis infection. This test should not be used as a primary diagnostic approach for syphilis. To screen for undiagnosed syphilis infection, order SYPH1 / Syphilis IgG with Reflex, Enzyme Immunoassay, Serum.

**Specimen Requirements:**  
**Supplies:** Sarstedt Aliquot Tube, 5 mL (T914)  
**Collection Container/Tube:**  
**Preferred:** Serum gel  
**Acceptable:** Red top  
**Submission Container/Tube:** Plastic vial  
**Specimen Volume:** 0.5 mL  
**Collection Instructions:** Centrifuge and aliquot serum into a plastic vial.

**Forms:**  
If not ordering electronically, complete, print, and send [Infectious Disease Serology Test Request](#) (T916) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Frozen	14 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
616863	RPR Screen w/ Reflex to Titer,S	Alphanumeric		20507-0

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**  
No

**CPT Code Information:**  
86592  
86593-Rapid Plasma Reagin Titer (if appropriate)

**Reflex Tests:**

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
RPRT2	RPR Titer,S	1	86593	No	No

# Test Definition: RPRT1

Rapid Plasma Reagin (RPR) Screen with Reflex  
to Titer, Serum

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
RPRT2	616865	RPR Titer,S	Alphanumeric		31147-2

Reference Values:

Negative  
Reference values apply to all ages.