

**Reporting Title:** HE4, S**Performing Location:** Rochester**Specimen Requirements:****Patient Preparation:** For 12 hours before specimen collection, patient **should not** take multivitamins or dietary supplements (eg, hair, skin, and nail supplements) containing biotin (vitamin B7).**Supplies:** Sarstedt Aliquot Tube, 5 mL (T914)**Collection Container/Tube:****Preferred:** Serum gel**Acceptable:** Red top**Submission Container:** Plastic vial**Specimen Volume:** 1 mL**Collection Instructions:** Centrifuge and aliquot serum into a plastic vial.**Forms:**If not ordering electronically, complete, print, and send an [Oncology Test Request](#) (T729) with the specimen.

| Specimen Type | Temperature        | Time    | Special Container |
|---------------|--------------------|---------|-------------------|
| Serum         | Refrigerated       | 7 days  |                   |
|               | Frozen (preferred) | 84 days |                   |

**Result Codes:**

| Result ID | Reporting Name | Type    | Unit   | LOINC®  |
|-----------|----------------|---------|--------|---------|
| HE4       | HE4, S         | Numeric | pmol/L | 55180-4 |

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

86305

**Reference Values:**

Females: &lt; or =140 pmol/L

Males: Not applicable