

**Reporting Title:** Thyroperoxidase Ab, S  
**Performing Location:** Rochester

**Specimen Requirements:**

**Patient Preparation:** For 12 hours before specimen collection do not take multivitamins or dietary supplements containing biotin (vitamin B7), which is commonly found in hair, skin, and nail supplements and multivitamins.

**Supplies:** Sarstedt 5 mL Aliquot Tube (T914)

**Collection Container/Tube:** Red top

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 0.6 mL

**Collection Instructions:** Centrifuge and aliquot serum into plastic vial within 2 hours of collection.

**Forms:**

Specimen Type	Temperature	Time	Special Container
Serum Red	Ambient	7 days	
	Refrigerated (preferred)	7 days	
	Frozen	30 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
TPO	Thyroperoxidase Ab, S	Numeric	IU/mL	8099-4

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

86376

**Reference Values:**

<9.0 IU/mL

Reference values apply to all ages.