

The following applies to IMMAU / Inborn Errors of Immunity with Immune Dysregulation and Autoimmunity Gene Panel, Varies. Testing is performed to evaluate for the presence of variants in coding regions and extending to +/- 10 base pairs of adjacent intronic sequence on either side of the coding exons of the genes analyzed. In addition, the analysis will cover select non-coding variants. Next-generation sequencing and/or a polymerase chain reaction-based quantitative method is performed to test for the presence of copy number variants in the genes analyzed. Confirmation of select reportable variants may be performed by alternate methodologies based on internal laboratory criteria.

This list is current from March 2023 to the present. This document is intended to highlight additional evaluations for variants of high clinical interest to providers as well as technical limitations. However, this document does not comprehensively reflect all genomic regions covered by this test. For questions regarding transcripts, genes or regions covered, contact the laboratory at 800-533-1710.

Genomic Build: GRCh37 (hg19) unless otherwise specified

| Gene          | Reference Transcript | Additional Evaluations   | Technical Limitations                          |
|---------------|----------------------|--|--|
| <i>AIRE</i>   | NM_000383.4          | -  | -  |
| <i>BACH2</i>  | NM_021813.4          | -  | -  |
| <i>CARD11</i> | NM_032415.6          | -  | -  |
| <i>CASP10</i> | NM_032977.4          | -  | -  |
| <i>CASP8</i>  | NM_001228.4          | -  | -  |
| <i>CD3G</i>   | NM_000073.2          | -  | -  |
| <i>CTLA4</i>  | NM_005214.5          | -  | -  |
| <i>DEF6</i>   | NM_022047.4          | -  | -  |
| <i>FADD</i>   | NM_003824.3          | -  | -  |
| <i>FASLG</i>  | NM_000639.3          | c.-261T>C  | -  |
| <i>FERMT1</i> | NM_017671.4          | c.-20A>G; c.1139+740G>A  | CNV analysis in exon 9 is not performed.       |
| <i>FOXP3</i>  | NM_014009.4          | c.-105_-23+1303del1386 (also known as g.del-6247_-4859); c.-7G>T; c.*876 to c.*878 | -  |
| <i>IL10</i>   | NM_000572.3          | -  | -  |
| <i>IL10RA</i> | NM_001558.3          | -  | -  |
| <i>IL10RB</i> | NM_000628.5          | -  | -  |
| <i>IL2RA</i>  | NM_000417.3          | -  | -  |
| <i>IL2RB</i>  | NM_000878.5          | -  | -  |
| <i>ITCH</i>   | NM_031483.7          | -  | -  |
| <i>JAK1</i>   | NM_002227.4          | -  | -  |
| <i>LRBA</i>   | NM_006726.4          | -  | CNV analysis in exons 20, 39 is not performed. |
| <i>ORAI1</i>  | NM_032790.3          | -  | -  |
| <i>PEPD</i>   | NM_000285.4          | -  | -  |
| <i>PRKCD</i>  | NM_006254.4          | -  | -  |
| <i>RIPK1</i>  | NM_003804.6          | c.688_688+20del  | -  |
| <i>STAT3</i>  | NM_139276.2          | c.1282-89C>T   | -  |
| <i>STAT5B</i> | NM_012448.4          | -  | CNV analysis in exons 6-8 is not performed.    |
| <i>STIM1</i>  | NM_003156.3          | -  | -  |
| <i>TET2</i>   | NM_001127208.2       | -  | -  |

*Targeted Genes and Methodology Details  
for Inborn Errors of Immunity with Immune  
Dysregulation and Autoimmunity Gene Panel* (continued)

| <b>Gene</b>  | <b>Reference Transcript</b> | <b>Additional Evaluations</b> | <b>Technical Limitations</b>              |
|--------------|-----------------------------|-------------------------------|---|
| <i>TGFB1</i> | NM_000660.7                 | -                             | -   |
| <i>TPP2</i>  | NM_003291.4                 | -                             | CNV analysis in exon 28 is not performed. |