

## Overview

### Useful For

Processing the bone marrow or peripheral blood specimen but delaying fluorescence in situ hybridization analysis while preliminary morphologic assessment is in process

### Reflex Tests

Test Id	Reporting Name	Available Separately	Always Performed
PHLDB	Probe, Each Additional (PHLDF)	No, (Bill Only)	No
BALPB	Probe, Each Additional (BALPF)	No, (Bill Only)	No
EOSMB	Probe, Each Additional (EOSMF)	No, (Bill Only)	No
EOS3B	Probe, Tri-color (EOSMF)	No, (Bill Only)	No
BALMB	Probe, Each Additional (BALMF)	No, (Bill Only)	No
BAL3B	Probe, Tri-color (BAL)	No, (Bill Only)	No
BLPMB	Probe, Each Additional (BLPMF)	No, (Bill Only)	No
BALAB	Probe, Each Additional (BALAF)	No, (Bill Only)	No
MFCDB	Probe, Each Additional (MFCDF)	No, (Bill Only)	No
EOSDB	Probe, Each Additional (EOSDF)	No, (Bill Only)	No
TLPDB	Probe, Each Additional (TLPDF)	No, (Bill Only)	No
COGBB	Probe, Each Additional (COGBF)	No, (Bill Only)	No
COGTB	Probe, Each Additional (COGTF)	No, (Bill Only)	No
COGMB	Probe, Each Additional (COGMF)	No, (Bill Only)	No
HEMMB	Probe, Each Additional (HEMMF)	No, (Bill Only)	No
MDSDB	Probe, Each Additional (MDSDF)	No, (Bill Only)	No
MDSMB	Probe, Each Additional (MDSMF)	No, (Bill Only)	No

AMLAF	Adult AML, FISH	Yes	No
AMLMF	AML, Specified FISH	Yes	No
AMLPF	Pediatric AML, FISH	Yes	No
BALAF	Adult ALL (B-cell), FISH	Yes	No
BALMF	ALL (B-cell), Specified FISH	Yes	No
BALPF	Pediatric ALL (B-cell), FISH	Yes	No
BLPMF	B-cell Lymphoma, Specified FISH	Yes	No
CILDF	Cong Infantile Leukemia, Diag FISH	Yes	No
CILMF	Cong Infantile Leukemia, Spec FISH	Yes	No
CLLDF	CLL, Diagnostic FISH	Yes	No
CLLMF	CLL, Specified FISH	Yes	No
COGBF	COG, ALL (B-cell), FISH	Yes	No
COGMF	COG, AML, FISH	Yes	No
COGTF	COG, ALL (T-cell), FISH	Yes	No
EOSDF	Chronic Eosinophilia, Diag FISH	Yes	No
EOSMF	Chronic Eosinophilia, Spec FISH	Yes	No
HEMMF	Hematologic Specified FISH	Yes	No
MDSDF	MDS, Diagnostic FISH	Yes	No
MDSMF	MDS, Specified FISH	Yes	No
MFCDF	Myeloma Fixed Cell, High Risk, FISH	Yes	No
PHLDF	Ph-like ALL(B-cell), Diag FISH	Yes	No
TALAF	Adult ALL (T-cell), FISH	Yes	No
TALMF	ALL (T-cell), Specified FISH	Yes	No
TALPF	Pediatric ALL (T-cell), FISH	Yes	No
TLPDF	T-cell Lymphoma B/BM, Diag FISH	Yes	No
TLPMF	T-cell Lymphoma B/BM, Spec FISH	Yes	No
TALAB	Probe, Each Additional (TALAF)	No, (Bill Only)	No
TALPB	Probe, Each Additional (TALPF)	No, (Bill Only)	No
TLPMB	Probe, Each Additional (TLPMF)	No	No
TALMB	Probe, Each Additional (TALMF)	No, (Bill Only)	No

AMLAB	Probe, Each Additional (AMLAF)	No, (Bill Only)	No
AMLMB	Probe, Each Additional (AMLMF)	No, (Bill Only)	No
AMLPB	Probe, Each Additional (AMLPF)	No, (Bill Only)	No
CLLDB	Probe, Each Additional (CLLDF)	No, (Bill Only)	No
CLLMB	Probe, Each Additional (CLLMF)	No, (Bill Only)	No

### Testing Algorithm

This test is designed to hold the sample and delay fluorescence in situ hybridization (FISH) testing while preliminary morphologic assessment or flow cytometry testing is in process.

**Hold policy:** Upon sample receipt, the specimen will be held in the laboratory. FISH testing will not be performed unless the client contacts the laboratory and indicates that FISH testing is desired. The client **must** contact the Cytogenetics Laboratory at 800-533-1710 by 4 p.m. (Central time) no later than 4 business days (96 hours) after the specimen was collected. If no notification is received by this time, the order will be processed as "canceled." Weekend communication can be deferred until Monday.

### Method Name

Direct Preparation of Specimen

### NY State Available

Yes

## Specimen

### Specimen Type

Varies

### Ordering Guidance

**This test is designed to hold blood and bone marrow specimens only.**

This test does not apply to chromosome analysis. If specimen is to be held for chromosome analysis, order HOLDC / Hematologic Disorders, Chromosome Hold, Varies.

Due to stability issues, test PCPDS / Plasma Cell Proliferative Disorder, High-Risk with Reflex Probes, Diagnostic FISH Evaluation, Bone Marrow **cannot be** added onto specimens held under this test.

### Shipping Instructions

Advise Express Mail or equivalent if not on courier service.

**Necessary Information**

Provide a reason for testing with each specimen and bone marrow pathology report (if available). The laboratory will not reject testing if this information is not provided, but appropriate testing and interpretation may be compromised or delayed.

**Specimen Required**

**Submit only 1 of the following specimens:**

**Preferred:**

**Specimen Type:** Bone marrow

**Container/Tube:**

**Preferred:** Yellow top (ACD)

**Acceptable:** Green top (sodium heparin), lavender top (EDTA)

**Specimen Volume:** 1 to 2 mL

**Collection Instructions:** Invert several times to mix bone marrow.

**Acceptable:**

**Specimen Type:** Blood

**Container/Tube:**

**Preferred:** Yellow top (ACD)

**Acceptable:** Green top (sodium heparin), lavender top (EDTA)

**Specimen Volume:** 6 mL

**Collection Instructions:** Invert several times to mix blood.

**Forms**

If not ordering electronically, complete, print, and send a [Hematopathology/Cytogenetics Test Request](#) (T726) with the specimen.

**Specimen Minimum Volume**

Blood: 2 mL

Bone marrow: 1 mL

**Reject Due To**

No specimen should be rejected.

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)		
	Standard		

**Clinical & Interpretive**

---

**Clinical Information**

Fluorescence in situ hybridization (FISH) analysis using gene-specific probes is a useful methodology to detect common, recurrent chromosome abnormalities for most hematologic malignancies. Based on morphologic review of the bone marrow or peripheral blood specimen by a hematopathologist, a determination of additional appropriate testing can be made. If the specimen does not show evidence of malignancy, FISH analysis may not be necessary. Depending on the diagnosis, conventional chromosome analysis may also be more informative.

**Reference Values**

Not applicable

**Interpretation**

If notified by the client, this test may be canceled, and a processing fee will be assessed.

If no notification to proceed with testing is received, this test will be reported as "cancelled."

**Cautions**

No significant cautionary statements

**Performance****Method Description**

The specimen will be held in the laboratory while preliminary morphologic assessment is in process to determine if fluorescence in situ hybridization (FISH) testing is appropriate.

**PDF Report**

No

**Day(s) Performed**

Monday through Sunday

**Report Available**

4 days

**Specimen Retention Time**

4 weeks

**Performing Laboratory Location**

Rochester

**Fees & Codes**

**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

Not Applicable

**CPT Code Information**

See individual reflex tests

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
HOLDF	Heme FISH Hold, B/BM	No LOINC Needed

Result ID	Test Result Name	Result LOINC® Value
51836	Result Summary	50397-9
51838	Interpretation	69965-2
CG666	Reason for Referral	42349-1
CG667	Specimen	31208-2
51839	Source	31208-2
CG791	Requested FISH Test	48767-8
51841	Method	85069-3
53433	Additional Information	48767-8
51842	Released by	18771-6