

Egg White IgG

Overview

Method Name

Enzyme Immunoassay (FEIA)

NY State Available

No

Specimen

Specimen Type Serum

Specimen Required

Collection Container/Tube:

Preferred: Red Top

Acceptable: Serum gel

Specimen Volume: 0.5 mL

Collection Instructions:

- 1. Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable.
- 2. Centrifuge and aliquot 0.5 mL of serum in a plastic vial.
- 3. Send refrigerate.

Specimen Minimum Volume

See Specimen Required

Reject Due To

Thawing: Cold OK; Warm OK	
---------------------------	--

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Ambient	7 days	
	Refrigerated (preferred)	28 days	
	Frozen	365 days	

Clinical & Interpretive

Reference Values



Egg White IgG

Reference ranges have not been established for food-specific IgG tests.

Interpretation

The clinical utility of food-specific IgG tests has not been established. These tests can be used in special clinical situations to select foods for evaluation by diet elimination and challenge in patients who have food-related complaints. It should be recognized that the presence of food-specific IgG alone cannot be taken as evidence of food allergy and only indicates immunologic sensitization by the food allergen in question. This test should only be ordered by physicians who recognize the limitations of the test.

Performance

PDF Report No

Day(s) Performed Monday through Friday

Report Available

5 to 7 days

Performing Laboratory Location

Eurofins Viracor

Fees & Codes

Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

86001

LOINC[®] Information

Test ID	Test Order Name	Order LOINC [®] Value
FEGWH	Egg White IgG	35535-4





Egg White IgG

Result ID	Test Result Name	Result LOINC [®] Value
FEGWH	Egg White IgG	35535-4