

**Overview****Method Name**

ImmunoCAP FEIA

**NY State Available**

Yes

**Specimen****Specimen Type**

Serum

**Specimen Required**

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

|           |                       |
|-----------|-----------------------|
| Hemolysis | Mild OK; Gross reject |
| Lipemia   | Mild OK; Gross reject |
| Icterus   | Mild OK; Gross reject |
| Other     | NA                    |

**Specimen Stability Information**

| Specimen Type | Temperature              | Time     | Special Container |
|---------------|--------------------------|----------|-------------------|
| Serum         | Ambient                  | 28 days  |                   |
|               | Refrigerated (preferred) | 28 days  |                   |
|               | Frozen                   | 365 days |                   |

**Clinical & Interpretive****Reference Values**

&lt;0.35 kU/L

**Interpretation**

| <u>Class</u> | <u>IgE (kU/L)</u> | <u>Comment</u>       |
|--------------|-------------------|----------------------|
| 0            | <0.10             | Negative             |
| 0/1          | 0.10 – 0.34       | Equivocal/Borderline |
| 1            | 0.35 – 0.69       | Low Positive         |
| 2            | 0.70 – 3.49       | Moderate Positive    |
| 3            | 3.50 – 17.49      | High Positive        |
| 4            | 17.50 – 49.99     | Very High Positive   |
| 5            | 50.00 – 99.99     | Very High Positive   |
| 6            | >99.99            | Very High Positive   |

## Performance

### PDF Report

No

### Day(s) Performed

Monday through Friday

### Report Available

1 to 6 days

### Performing Laboratory Location

Eurofins Viracor

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### CPT Code Information

86003

### LOINC® Information

| Test ID | Test Order Name     | Order LOINC® Value |
|---------|---------------------|--------------------|
| FNBE    | Bean Navy/White IgE | 6279-4             |

| Result ID | Test Result Name    | Result LOINC® Value |
|-----------|---------------------|---------------------|
| Z3992     | Bean Navy/White IgE | 6279-4              |

---

|       |       |          |
|-------|-------|----------|
| Z3993 | CLASS | 102173-2 |
|-------|-------|----------|