

# **Test Definition: FPOLO**

Poliovirus (Types 1, 3) Antibodies, Neutralization

#### Overview

#### **Method Name**

Culture/Neutralization

#### **NY State Available**

Yes

## **Specimen**

## **Specimen Type**

Serum

## **Specimen Required**

**Specimen Type:** Serum **Container/Tube:** Red or SST **Specimen Volume:** 1 mL

**Collection Instructions:** Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum refrigerate in a plastic vial.

## **Specimen Minimum Volume**

0.5 mL

#### **Reject Due To**

## **Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Ambient	5 days	
	Refrigerated (preferred)	14 days	
	Frozen	30 days	

## **Clinical & Interpretive**

#### **Clinical Information**

This sensitive procedure is recommended for vaccine response testing and type-specific serodiagnosis of recent poliovirus infection. It can also serve as an aid for diagnosing immune deficiency disorders.



## **Test Definition: FPOLO**

Poliovirus (Types 1, 3) Antibodies, Neutralization

#### **Reference Values**

Polio 1 Titer: <1:8 Polio 3 Titer: <1:8

The presence of neutralizing serum antibodies (titers 1:8 up to >1:128) against polioviruses implies lifelong immunity. Some persons without detectable titers (<1:8) may also be immune as demonstrated by elicitation of a secondary-type serum antibody response upon rechallenge with live polio vaccine.

#### **Performance**

### **PDF Report**

No

## Day(s) Performed

Monday and Thursday

#### Report Available

1 to 10 days

## **Performing Laboratory Location**

**Quest Diagnostics** 

#### Fees & Codes

#### Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

#### **Test Classification**

This test was developed and its analytical performance characteristics have been determined by Quest Diagnostics. It has not been cleared or approved by FDA. This assay has been validated pursuant to the CLIA regulations and is used for clinical purposes.

#### **CPT Code Information**

86382 x 2

#### **LOINC®** Information

|--|



# **Test Definition: FPOLO**

Poliovirus (Types 1, 3) Antibodies, Neutralization

FPOLO	Poliovirus (Types 1, 3) Ab, Neut	68320-1
-------	----------------------------------	---------

Result ID	Test Result Name	Result LOINC® Value
Z4806	Polio 1 Titer	22446-9
Z4807	Polio 3 Titer	22450-1