

Overview**Reflex Tests**

| Test Id | Reporting Name | Available Separately | Always Performed |
|---------|------------------------|----------------------|------------------|
| FCOKR | Cocaine, MS, WB/SP Rfx | No | No |

Testing Algorithm

If the Cocaine Analysis, whole blood is positive, then the Cocaine, Confirmation (FCOKR) will be performed at an additional charge.

Method Name

Immunoassay (IA)

Gas Chromatography with mass spectrometry (GC/MS); confirmation

NY State Available

Yes

Specimen**Specimen Type**

Whole blood

Specimen Required

Draw blood in a purple-top (EDTA), green-top (sodium-heparin) or grey-top (NaF/oxalate) tube(s). Send 7 mL of whole blood refrigerated in a plastic vial.

Specimen Minimum Volume

3 mL

Reject Due To

| | |
|-----------|----|
| Hemolysis | NA |
| Lipemia | NA |
| Icterus | NA |
| Other | NA |

Specimen Stability Information

| Specimen Type | Temperature | Time | Special Container |
|---------------|-------------|------|-------------------|
|---------------|-------------|------|-------------------|

| | | | |
|-------------|--------------------------|----------|--|
| Whole blood | Refrigerated (preferred) | 14 days | |
| | Frozen | 180 days | |
| | Ambient | 72 hours | |

Clinical & Interpretive

Reference Values

Cocaine Screen:

Cocaine and Metabolite, UA – Negative; Cutoff: 25 ng/mL

Cocaine Confirmation:

Cocaine

Benzoylcegonine

Confirmation threshold: 10 ng/mL

Performance

PDF Report

No

Day(s) Performed

Monday through Sunday

Report Available

5 to 9 days

Performing Laboratory Location

Medtox Laboratories, Inc.

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.

CPT Code Information

80307

80353-if applicable

LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---------|-------------------------|--------------------|
| FCOKE | Cocaine, Scr w/Conf, WB | 16633-0 |

| Result ID | Test Result Name | Result LOINC® Value |
|-----------|------------------------|---------------------|
| Z4826 | Cocaine and Metabolite | 73823-7 |