

## Overview

**Method Name**

ImmunoCAP FEIA

**NY State Available**

Yes

## Specimen

**Specimen Type**

Serum

**Specimen Required****Specimen Type:** Serum**Container/Tube:** Red or SST**Specimen Volume:** 0.5 mL**Collection Instructions:** Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

Thawing:	Warm OK; Cold OK
Gross lipemia:	Reject

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

## Clinical & Interpretive

**Reference Values**

---

<0.35 kU/L

**Interpretation**

Class	IgE (kU/L)	Comment
0	<0.10	Negative
0/1	0.10-0.34	Equivocal/Borderline
1	0.35 - 0.69	Low Positive
2	0.70 - 3.49	Moderate Positive
3	3.50 - 17.49	High Positive
4	17.50 - 49.99	Very High Positive
5	50.00 - 99.99	Very High Positive
6	>99.99	Very High Positive

**Performance****PDF Report**

No

**Day(s) Performed**

Monday through Friday

**Report Available**

3 to 6 days

**Performing Laboratory Location**

Eurofins Viracor

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

**CPT Code Information**

---

86003

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
FFTMV	Trichophyton ment var interdig IgE	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z5619	Trichophyton ment var interdig IgE	25801-2
Z5620	Class	102906-5