

Overview

Method Name

Gas Chromatography (GC)

NY State Available

Yes

Specimen

Specimen Type

Varies

Specimen Required

*****Must submit one specimen per order. Specimens cannot be shared between multiple orders.*****

Submit only 1 of the following specimens:

Serum

Draw blood in a plain, red-top tube(s). **Serum gel tube is not acceptable.** Centrifuge and send 2 mL of serum refrigerate in plastic vial.

Note: 1. Indicate serum on request form.

2. Label specimen appropriately (serum).

Plasma

Draw blood in a lavender-top (EDTA) tube(s). **Plasma gel tube is not acceptable.** Centrifuge and send 2 mL of EDTA plasma refrigerate in plastic vial.

Note: 1. Indicate plasma on request form.

2. Label specimen appropriately (plasma).

Specimen Minimum Volume

1.0 mL

Reject Due To

Hemolysis	NA
Lipemia	NA
Icteric	NA
Other	Polymer gel separation tube (SST or PST)

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated		

Clinical & Interpretive

Reference Values

Reporting limit determined each analysis

Synonym(s): Diprivan

Patients required a mean blood propofol concentration of 4.05 +/- 1.01 mcg/mL for major surgery and 2.97 +/- 1.07 mcg/mL for non-major surgery. Blood propofol concentrations at which 50% of patients were awake and oriented after surgery were 1.07 and 0.95 mcg/mL respectively.

Psychomotor performance returned to baseline at blood propofol concentrations of 0.38 - 0.43 mcg/mL.

Performance

PDF Report

No

Day(s) Performed

Monday through Sunday

Report Available

8 to 12 days

Specimen Retention Time

2 weeks

Performing Laboratory Location

NMS Labs

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

80375

G0480 (if appropriate)

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FPROP	Propofol	34673-4

Result ID	Test Result Name	Result LOINC® Value
Z1448	Propofol	34673-4
Z1869	Reporting Limit	19147-8