

TEST OBSOLETE

Notification Date: August 30, 2023 Effective Date: August 31, 2023

17-Hydroxyprogesterone

Test ID: F17HY

Explanation:

Effective August 31, 2023, Test ID F17HY, referred to Quest Diagnostics SJC Laboratory, will be obsolete due to availability of Mayo Test ID OHPG.

Recommended Alternative Test: 17-Hydroxyprogesterone, Serum

Test ID: OHPG

Methodology:

Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS)

Reference Values:

Children: Preterm infants Preterm infants may exceed 630 ng/dL, however, it is uncommon to see levels reach 1,000 ng/dL. Term infants 0-28 days: <630 ng/dL Levels fall from newborn (<630 ng/dL) to prepubertal gradually within 6 months. Prepubertal males: <110 ng/dL Prepubertal females: <100 ng/dL

Adults: Males: <220 ng/dL Females Follicular: <80 ng/dL Luteal: <285 ng/dL Postmenopausal: <51 ng/dL

Note: For pregnancy reference ranges, see: Soldin OP, Guo T, Weiderpass E, et al: Steroid hormone levels in pregnancy and 1 year postpartum using isotope dilution tandem mass spectrometry. Fertil Steril. 2005 Sept;84(3):701-710

Specimen Requirements:

| Container/Tube: | Red Top |
|---------------------------------|--|
| Specimen Volume: | 0.6 mL |
| Collection Instructions: | Centrifuge and aliquot serum into a plastic vial |
| Minimum Volume: | 0.25mL |

Specimen Stability Information:

| Specimen Type | Temperature | Time |
|---------------|--------------------------|---------|
| Serum Red | Refrigerated (preferred) | 28 days |
| | Ambient | 7 days |
| | Frozen | 28 days |

CPT Code: 83498

Day(s) Setup: Monday through Friday

Report Available: 2 to 5 days

Questions

Contact MCL Referrals Supervisor Amy Bluhm at 800-533-1710.