MAYO CLINIC LABORATORIES

TEST STATUS - TEST DOWN

Notification Date: December 14, 2023 Effective Date: Immediately

Methadone and Metabolites, Serum

Test ID: MDNS

Explanation: This test will be made non-orderable effective immediately due to instrumentation issues and assay redevelopment. An update will be provided when testing resumes.

Recommended Alternative Test:

Methadone and Metabolite, Serum

Test ID: FMTMS

Methodology: High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)

Reference Values:

Reporting limits(s) determined each analysis.

Methadone: None Detected ng/mL EDDP: None Detected ng/mL

Specimen Requirements:

Must submit one specimen per order. Specimens cannot be shared between multiple orders.

Specimen Type: Serum

Collection Container/Tube: Red top

Submission Container/Tube: Plastic, preservative-free vial

Specimen Volume: 1 mL Collection Instructions:

1. Draw blood in a red top tube(s). Serum gel tube is not acceptable.

2. Centrifuge and send 1 mL of serum refrigerated in a plastic, preservative-free vial.

Specimen Stability Information:

| Specimen Type | Temperature | Time |
|---------------|-------------------------|----------|
| Serum Red | Refrigerate (preferred) | 29 days |
| | Frozen | 365 days |

CPT Code:

80358

G0480 (if appropriate)

Miscellaneous Referral Test Fee: \$158.00

Day(s) Performed: Monday through Sunday Report Available: 6 to 10 days

Questions

Contact Brandon DeBoom, Laboratory Resource Coordinator at 800-533-1710.