



Allo-isolucine, Blood Spot

Test ID: ALLOI

Explanation: On the effective date the Specimen Stability Time will be updated.

Current Specimen Stability			
Specimen	Temperature	Time	
Whole blood	Ambient (preferred)		FILTER PAPER
	Frozen		FILTER PAPER
	Refrigerated		FILTER PAPER

New Specimen Stability			
Specimen	Temperature	Time	Special Container
Whole blood	Ambient (preferred)	59 days	FILTER PAPER
	Frozen	59 days	FILTER PAPER
	Refrigerated	59 days	FILTER PAPER

Questions

Contact Melissa Tricker-Klar, Laboratory Technologist Resource Coordinator at 800-533-1710.