



Amylase, Isoenzymes, Serum

Test ID: AMISO

Explanation: On the effective date, the preferred temperature will change from Frozen to Refrigerated.

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen (preferred)	30 days
	Ambient	7 days
	Refrigerated	30 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	30 days
	Ambient	7 days
	Frozen	30 days

Questions

Contact Nancy Benson, Laboratory Resource Coordinator at 800-533-1710.