



---

## Sequential Maternal Screening, Part 1, Serum

Test ID: SEQA

### Explanation:

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	14 days
	Ambient	7 days
	Frozen	90 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	7 days
	Ambient	7 days
	Frozen	90 days

### Questions

Contact Bethany Feind, Laboratory Resource Coordinator at 800-533-1710.