

## Tyrosinemia Follow Up Panel, Self-Collect, Blood Spot

**Test ID:** TYRSC

**Explanation:** The assay performance issues affecting this test have been resolved. Testing will resume on the effective date.

Specimen Stability, Reference Values and CPT codes have been updated, as indicated below.

Current Specimen Stability			
Specimen	Temperature	Time	Special Container
Whole blood	Ambient (preferred)	90 days	FILTER PAPER
	Frozen	90 days	FILTER PAPER
	Refrigerate	90 days	FILTER PAPER

New Specimen Stability			
Specimen	Temperature	Time	Special Container
Whole blood	Ambient (preferred)	21 days	FILTER PAPER
	Frozen	10 days	FILTER PAPER
	Refrigerate	10 days	FILTER PAPER

Current Reference Value
TYROSINE: <4 weeks 40.0-280.0 nmol/mL > or =4 weeks 25.0-150.0 nmol/mL  PHENYLALANINE: 27.0-107.0 nmol/mL  METHIONINE: 11.0-45.0 nmol/mL  SUCCINYLACETONE: <1.00 nmol/mL  NITISINONE: <0.7 nmol/mL

New Reference Value
TYROSINE: <4 weeks 40.0-280.0 nmol/mL > or =4 weeks 25.0-150.0 nmol/mL  PHENYLALANINE: 27.0-107.0 nmol/mL  SUCCINYLACETONE: <b>&lt;1.6 nmol/mL</b>  NITISINONE: <b>&lt;0.6 nmol/mL</b>

<b>Current CPT Code</b>
84510
84030
82131
82542
80299
82542 only (if appropriate for government payers)

<b>New CPT Code</b>
84510
84030
82542
80299
82542 only (if appropriate for government payers)

## **Questions**

Contact Melissa Tricker-Klar, Laboratory Resource Coordinator at 800-533-1710.