

# Ma2 Antibody by ELISA, Spinal Fluid

Test ID: MA2EC

## Useful for:

Evaluating patients with suspected paraneoplastic encephalitides using spinal fluid specimens

### Methods:

Enzyme-Linked Immunosorbent Assay (ELISA)

## **Reference Values:**

Negative

## **Specimen Requirements:**

Container/Tube:	Sterile vial	
Preferred:	Vial number 1	
Acceptable:	Any vial	
Specimen Volume:	2 mL	
Minimum Volume:	1 mL	

#### Note:

Provide the following information:

1.Relevant clinical information

2. Ordering provider name, phone number, mailing address, and e-mail address

# **Specimen Stability Information:**

Specimen Type	Temperature	Time
CSF	Refrigerated (preferred)	28 days
	Ambient	72 hours
	Frozen	28 days

# Cautions:

A negative Ma2 antibody test result does not exclude autoimmune neurological disease or cancer.

#### **CPT Code:**

83516

Day(s) Performed: Tuesday, Friday

Report Available: 3 to 5 days

## Note:

The following referral test code(s) will become obsolete.

Test Name	Test ID	Referral Lab Code	Referral Lab
Recombx MaTa Autoantibody Test	FRMTA	4683	Athena Diagnostics

# Questions

Contact Amy Ennis or Steven Monson, Laboratory Resource Coordinator at 800-533-1710.