



Des-Gamma-Carboxy Prothrombin, Serum

Test ID: DCP

Explanation: Additional stability added and preferred sample type changes are reflected below.

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen (preferred)	90 days
	Refrigerated	7 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen (preferred)	90 days
	Refrigerated	7 days
	Ambient	3 days

Current Specimen Required	
Preferred	Red Top
Acceptable	Serum Gel

New Specimen Required	
Preferred	Serum Gel
Acceptable	Red Top

Questions

Contact Bethany Feind, Laboratory Resource Coordinator at 800-533-1710.