

TEST CHANGE-UPDATE

Notification Date: April 28, 2023 Effective Date: immediately

Attention: Previously announced specimen stability, effective April 28, 2023, will not be implemented. Specimen stability will remain as below.

Entamoeba histolytica Antibody, Serum

Test ID: EHOLG

Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen (preferred)	30 days
	Refrigerated	7 days

Questions

Contact Dunisha Messmer, Laboratory Resource Coordinator at 800-533-1710.