

Thyroglobulin, Tumor Marker Reflex, Serum

Test ID: HTGR; performed at Mayo Clinic Laboratories Florida.

Explanation: on the effective date specimen volume and stability will be updated as indicated in the tables below.

Current Requested Volume
1 mL

New Requested Volume
2 mL

Current Specimen Minimum Volume
1 mL

New Specimen Minimum Volume
1.25 mL

Current Specimen Stability		
Specimen	Temperature	Time
Serum Red	Refrigerated (preferred)	7 days
	Frozen	30 days

New Specimen Stability		
Specimen	Temperature	Time
Serum Red	Refrigerated (preferred)	7 days
	Ambient	72 hours
	Frozen	30 days

Questions

Contact Bonnie Meyers, Laboratory Resource Coordinator at 800-533-1710.