

**Explanation:** We have been notified by Eurofins Viracor of the following reference range, and stability updates for the test ID's listed below which are effective July 24, 2023. New sample reports will be available after the effective date.

Test ID	Test Name	Current Reference Value	New Reference Value
FIL4S	Interleukin-4 (IL-4) Serum	<2.0	<b>&lt;0.5</b>
FIL8S	Interleukin-8 (IL-8) Serum	<57.8	<b>&lt;12.4</b>
FIFNY	Interferon-gamma (IFN-γ) Serum	<2.0	<b>&lt;8.6</b>
FNFAS	TNF-alpha (TNF-a) Serum	<5.6	<b>&lt;1.4</b>

**Test ID: FIL8S**

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen	28 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen	<b>14 days</b>

**Test IDs: FIL4S, and FNFAS**

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen	365 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen	<b>14 days</b>

**Questions:**

Contact MCL Referrals Supervisor Amy Bluhm at 800-533-1710.