

# Heartland Virus, RNA, Molecular Detection, PCR, Spinal Fluid

Test ID: HRTVC

### Useful for:

Aiding in the diagnosis of central nervous system infection caused by Heartland virus using spinal fluid specimens.

#### Methods:

Real-Time Polymerase Chain Reaction (PCR)

#### **Reference Values:**

Negative

#### **Specimen Requirements:**

Container/Tube: Sterile vial Specimen Volume: 1 mL Collection Instructions: 1. Send specimen from collection vial 2.

## 2. Do not centrifuge or heat inactivate.

#### **Specimen Stability Information:**

Specimen Type	Temperature	Time
CSF	Refrigerated (preferred)	14 days
	Frozen	14 days

#### Cautions:

Results should be used in conjunction with clinical presentation and exposure history.

Negative Heartland virus (HRTV) polymerase chain reaction results may occur if the specimen was collected more than 7 days following symptom onset. Serologic testing for the presence of antibodies to HRTV may be recommended in such cases.

**CPT Code:** 87798

Day(s) Performed: Monday through Friday

Report Available: Same day/1 to 5 days

# Questions

Contact James Conn, Laboratory Resource Coordinator at 800-533-1710.