



Procollagen I Intact N-Terminal, Serum

Test ID: PINP

Explanation:

Reagent backorder issues have been resolved. Testing for PINP: Procollagen I Intact N-Terminal, Serum will resume immediately.

Note:

Stability will be updated as noted below.

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen (preferred)	28 days
	Ambient	7 days
	Refrigerated	7 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen (preferred)	180 days
	Ambient	7 days
	Refrigerated	7 days

Questions

Contact Bethany Feind, Laboratory Resource Coordinator at 800-533-1710.