

Fungitell, CSF

Test ID: FUNSF

Explanation: Per review of test information for test FUNSF performed by Eurofins Viracor, the following change will be made effective immediately.

Current Specimen Stability		
Specimen	Temperature	Time
CSF	Frozen (preferred)	
	Ambient	4 days
	Refrigerated	7 days

New Specimen Stability		
Specimen	Temperature	Time
CSF	Frozen (preferred)	30 days
	Ambient	4 days
	Refrigerated	7 days

Questions:

Contact MCL Referrals Supervisor Amy Bluhm at 800-533-1710.