

Mycoplasma (Mycoplasmoides) pneumoniae with Macrolide Resistance Reflex, Molecular Detection, PCR, Varies

Test ID: MPRP

Explanation: Respiratory fluid and body fluid specimens submitted for this test will no longer be accepted if sent in viral transport media (VTM). These fluids should be transported in a sterile container not containing media.

| Current Fluid Specimen Requirement | New Fluid Specimen Requirement |
|---|--|
| <p>Specimen Type: Respiratory Supplies: M4-RT (T605) Sources: Bronchial washing, bronchoalveolar lavage, tracheal secretions, sputum Container/Tube: Preferred: Sterile container Acceptable: Specimen in M4, M4-RT, M5, M6, or universal transport medium Specimen Volume: 1 mL</p> <p>Specimen Type: Fluid Sources: Pleural, pericardial, cerebrospinal Container/Tube: Sterile vial Specimen Volume: 0.5 mL</p> | <p>Specimen Type: Respiratory Sources: Bronchial washing, bronchoalveolar lavage, tracheal secretions, sputum Container/Tube: Sterile container Specimen Volume: 1 mL</p> <p>Specimen Type: Fluid Sources: Pleural, pericardial, cerebrospinal Container/Tube: Sterile container Specimen Volume: 0.5 mL</p> |

Questions

Contact James Conn, Laboratory Resource Coordinator at 800-533-1710.