



Hemosiderin, Random, Urine

Test ID: UHSD2

Useful for:

Detecting hemosiderinuria, secondary to excess hemolysis, as in incompatible blood transfusions, severe acute hemolytic anemia, or hemochromatosis for external patients.

Method:

Rous Method

Reference Value:

Negative

Specimen Requirements:

- Container/Tube:** Plastic Urine Container
- Specimen Volume:** 13 mL
- Collection Instructions:**
1. Collect a random urine specimen
 2. No preservative
- Minimum Volume:** 12 mL

Specimen Stability Information:

Specimen Type	Temperature	Time
Urine	Refrigerated (preferred)	7 days
	Ambient	2 hours
	Frozen	7 days

CPT Code:

83070

Day(s) Performed: Monday through Sunday

Report Available: 1 day

Questions

Contact Nancy Benson, Laboratory Resource Coordinator at 800-533-1710.