



## Eculizumab, Serum

Test ID: ECULI

### Explanation:

Due to the implementation of new panel tests, the following specimen changes will occur on the effective date.

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (Preferred)	28 days
	Ambient	4 days
	Frozen	30 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen (Preferred)	30 days
	Ambient	4 days
	Refrigerated	28 days

### Questions

Contact Joshua Yang, Laboratory Resource Coordinator at 800-533-1710.