



Adenovirus DNA, Quantitative Real-Time PCR, Plasma

Test ID: FADNA

Explanation: Per notification from Quest Diagnostics, effective January 8, 2024, yellow-top (ACD) and PPT (white-top) tube types are no longer acceptable.

| Current Specimen Required |
|---|
| Plasma collected in yellow-top (ACD), Lavender-top (EDTA) or PPT (white-top) tube |

| New Specimen Required |
|---|
| Plasma collected in Lavender-top (EDTA) |

Questions:

Contact MCL Referrals Supervisor Amy Bluhm at 800-533-1710.