



Plasma Cell Proliferative Disorder, FISH, Tissue

Test ID: PLASF

Explanation:

On the effective date the requested volume for sending tissue slides will be updated as seen below. Please note that the preferred specimen type, tissue block, will remain unchanged.

| Current Requested Volume |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Acceptable: Tissue slides |
| Collection Instructions: For each probe set ordered, 2 consecutive, unstained, 5 micron-thick sections placed on positively charged slides, and 1 hematoxylin and eosin-stained slide. |

| New Requested Volume |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Acceptable: Tissue slides |
| Collection Instructions: Submit 10 consecutive unstained, positively charged, unbaked slides with 5 micron-thick sections of the tumor tissue and 1 slide stained with hematoxylin and eosin. |

Questions

Contact Joshua Couchene, Laboratory Resource Coordinator at 800-533-1710.