

HIV-1 and HIV-2 Antigen and Antibody Routine Screen, Plasma

Test ID: HVCOP

Explanation: Due to a method change, the following changes will be implemented on the effective date.

Current Method
Chemiluminescent Microparticle Immunoassay

New Method
Electrochemiluminescence Immunoassay (ECLIA)

Current Specimen Stability		
Specimen	Temperature	Time
Plasma	Refrigerated (preferred)	6 days
	Frozen	30 days

New Specimen Stability		
Specimen	Temperature	Time
Plasma EDTA	Frozen (preferred)	30 days
	Refrigerated	6 days

Current Report Available
1 to 2 days

New Report Available
1 to 3 days

Questions

Contact Dunisha Messmer, Laboratory Resource Coordinator at 800-533-1710.