

Hepatitis B e Antigen and Hepatitis B e Antibody, Serum

Test ID: HEAG

Explanation: Due to a method change, the following changes will be implemented on the effective date.

Current Method
Chemiluminescence Immunoassay

New Method
Electrochemiluminescence Immunoassay (ECLIA)

Current Specimen Stability		
Specimen	Temperature	Time
Serum SST	Frozen (preferred)	28 days
	Ambient	24 hours
	Refrigerated	7 days

New Specimen Stability		
Specimen	Temperature	Time
Serum SST	Frozen (preferred)	90 days
	Ambient	72 hours
	Refrigerated	6 days

Current Report Available
Same day/1 to 2 days

New Report Available
Same day/1 to 3 days

Questions

Contact Dunisha Messmer, Laboratory Resource Coordinator at 800-533-1710.