

Hepatitis A IgM Antibody, Serum

Test ID: HAIGM

Explanation: Due to a method change, the following changes will be implemented on the effective date.

Current Method	
Chemiluminescent Microparticle Immunoassay (CMIA)	

New Method	
Electrochemiluminescence Immunoassay (ECLIA)	

Current Requested Volume	
Volume	1 mL
Minimum Volume	0.4 mL

New Requested Volume	
Volume	0.6 mL
Minimum Volume	0.6 mL

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	7 days
	Ambient	72 hours

New Specimen Stability		
Specimen	Temperature	Time
Serum SST	Frozen (preferred)	90 days
	Ambient	72 hours
	Refrigerated	6 days

Questions

Contact Dunisha Messmer, Laboratory Resource Coordinator at 800-533-1710.