

Chromogranin A, Serum

Test ID: CGAK

Explanation: On the effective date, the following changes will be made to the Chromogranin A, Serum (CGAK) assay.

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen (preferred)	90 days
	Ambient	7 days
	Refrigerated	24 hours

New Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen (preferred)	90 days
	Ambient	48 hours
	Refrigerated	48 hours

Questions

Contact Bethany Feind, Laboratory Resource Coordinator at 800-533-1710.