

## Inhibin B, Serum

**Test ID:** INHB

**Explanation:** On the effective date, the following changes will be made to the Inhibin B, Serum (INHB) assay.

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	14 days
	Frozen	90 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	14 days
	Ambient	7 days
	Frozen	90 days

Current Requested Volume
1.0 mL

New Requested Volume
0.5 mL

### Questions

Contact Bethany Feind, Laboratory Resource Coordinator at 800-533-1710.