

## Inhibin A and B, Tumor Marker, Serum

**Test ID:** INHAB

**Explanation:** On the effective date, the following changes will be made to the Inhibin A and B, Tumor Marker, Serum (INHAB) assay.

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	14 days
	Frozen	90 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	14 days
	Ambient	7 days
	Frozen	90 days

### INHIBIN A, Tumor Marker

Current Reference Value	
Males	<2.0 pg/mL
Females	
<11 years	<4.7 pg/mL
11-17 years	<97.5 pg/mL
Premenopausal	<97.5 pg/mL
Postmenopausal	<2.1 pg/mL

New Reference Value	
Males	<5.0 pg/mL
Females	
<11 years	<5.0 pg/mL
11-17 years	<98 pg/mL
Premenopausal	<98 pg/mL
Postmenopausal	<5.0 pg/mL

Current Requested Volume
1.5 mL

New Requested Volume
1.1 mL

### Questions

Contact Bethany Feind, Laboratory Resource Coordinator at 800-533-1710.