

3-Methylglutaconic Aciduria Panel, Varies

Test ID: 3MGAP

Explanation: On the effective date, this assay will no longer accept *any anticoagulant* for the whole blood specimen type and will be limited to the anticoagulants shown below.

Current Specimen Required		New Specimen Required	
Specimen Type:	Whole blood	Specimen Type:	Whole blood
Preferred:	Lavender top (EDTA) or yellow top (ACD)	Required:	Lavender top (EDTA) or yellow top (ACD)
Acceptable:	Any anticoagulant	Specimen Volume:	3 mL
Specimen Volume:	3 mL	Collection Instructions:	1. Invert several times to mix blood 2. Send whole blood specimen in original tube. Do not aliquot.
Collection Instructions:	1. Invert several times to mix blood 2. Send whole blood specimen in original tube. Do not aliquot.	Specimen Stability Information:	Ambient (preferred) 4 days/Refrigerated 14 days
Specimen Stability Information:	Ambient (preferred) 4 days/Refrigerated 14 days	Minimum Volume:	1 mL
Minimum Volume:	1 mL		

Questions

Contact Michelle Rath, Laboratory Resource Coordinator at 800-533-1710.