

## Hepatitis B Virus Core Total Antibodies, with Reflex to Hepatitis B Virus Core IgM Antibody, Serum

**Test ID:** CORAB; performed at Mayo Clinic Laboratories Florida.

**Explanation:** Due to manufacturer's modification of stability times, the specimen stability information will be updated as indicated below.

Current Specimen Stability		
Specimen	Temperature	Time
Serum SST	Frozen (preferred)	90 days
	Ambient	72 hours
	Refrigerated	6 days

New Specimen Stability		
Specimen	Temperature	Time
Serum SST	Frozen (preferred)	90 days
	Refrigerated	6 days

### Questions

Contact Bonnie Meyers, Laboratory Resource Coordinator at 800-533-1710.