

Endothelin I

Explanation: The following new orderable, referred to Inter Science Institute, will be available immediately.

Test ID: FENDI

Methods:

Direct RIA

Reference Values:

4.0-9.0 pg/mL

Specimen Requirements:

Patient Preparation:

Patient should NOT be on any ACTH, Corticosteroids, or hypertension medications, if possible, for at least 48 hours prior to collection of specimen.

Submit only one of the following:

Container/Tube: Red or SST

Submission Container/Tube: 12x75 mm screw capped vial

Specimen Volume: 3 mL

Collection Instructions:

1. Draw blood in a plain, red-top or serum-gel tube(s).
2. Separate immediately
3. Send 3 mL of serum frozen in a plastic vial.

Minimum Volume: 1 mL

Or

Container/Tube: EDTA

Submission Container/Tube: 12x75 mm screw capped vial

Specimen Volume: 3 mL

Collection Instructions:

1. Draw blood in a lavender top (EDTA) tube(s).
2. Separate immediately
3. Send 3 mL of plasma frozen in a plastic vial.

Minimum Volume: 1 mL

Specimen Stability Information:

| Specimen Type | Temperature | Time |
|---------------|--------------------|----------|
| Varies | Frozen (preferred) | 180 days |
| | Refrigerated | 7 days |

Fee:
\$260.00

CPT Code:
83519

Day(s) Setup: Monday through Friday

Report Available: 12 to 16 days

Questions:

Contact MCL Referrals Supervisor Amy Bluhm at 800-533-1710.