
Spirochete Immunostain, Technical Component Only

Test ID: SPIRO

Explanation: Beginning on November 25th, 2024, in order to facilitate application of control tissue to the slide, test SPIRO will require a paraffin-embedded tissue block as the only acceptable specimen type. Unstained slides will no longer be accepted and testing will be canceled.

Current Specimen Required	New Specimen Required
<p>Supplies: Immunostain Technical Only Envelope (T693) Specimen Type: Tissue Container/Tube: Immunostain Technical Only Envelope Preferred: 2 Unstained positively charged glass slide (25- x 75- x 1-mm) per test ordered; sections 4-microns thick Acceptable: Formalin-fixed, paraffin-embedded (FFPE) tissue block</p>	<p>Supplies: Immunostain Technical Only Envelope (T693) Specimen Type: Tissue Container/Tube: Immunostain Technical Only Envelope Preferred: Formalin-fixed, paraffin-embedded (FFPE) tissue block</p>

Questions

Contact Richard Einerson/Sarah Turner, Laboratory Resource Coordinator at 800-533-1710.